

三通阀



Features:

- 1.Handle is stopped every 45 degrees angle
手柄45° 转角停顿
- 2.Luer lock rotates 360° to prevent tubing kinking
鲁尔螺旋管360° 无缠绕, 防止折管


Item	Color	Technical Data			Certificates
		Pressure resistance	Flow rate	Material	
SK45-R	Red	200KPa	Flow rate > 1000ml/10mins	PC& PE& POM& PP	CE certified
SK45-Y	Yellow				
SK45-B	Blue				



SK45-R

SK45-Y

SK45-B




SK46



SK47

成品输液套装—带三通阀



Extension Set

- Extension set for infusion and transfusion.
- Male luer lock
 - Female luer lock
 - 3 way stopcock
 - Tubing Length: 15/30/60/90/120cm



Gravity IV Set

- Standard IV set for infusion of drugs and liquids
- Drop chamber, 20 drops/ml, filter 15um
 - Roller clamp
 - Male luer lock
 - 3 way stopcock
 - Male luer lock



Pump IV Set

- Pump IV set for infusion of drugs and liquids
- Drop chamber, 20 drops/ml, filter 15um
 - Roller clamp
 - Male luer lock
 - 3 way stopcock
 - Needle free valve
 - Male luer lock



IV Burette Set

- Standard IV set for Burette set
- Vented spike
 - Roller clamp
 - 100/150ml Burette chamber(optional)
 - Roller clamp
 - Male luer lock
 - 3 way stopcock
 - Needle free valve
 - Male luer lock

