

成品输液器套装--止液型

Liquid Lock

- Prevent air infusion
- Higher filtration
- Efficiency of 8µm

Hydrophilic PES Membrane

Roller Clamp

- Average adjustment range: 19mm
- Hold the tube end
- Protect from sharp spike

Vented air cap

- Hydrophobic membrane 0.5µm

Drip chamber

- 20 drops/ml
- Φ16.3*72mm
- PVC/TPE

Prime Stop Cap

- Air pass
- No drops

Tubing

- DEHP free
- PVC free

A7F

64空气塞无滤网滴瓶头部, 20滴
64mm, Vented assembled chamber, 20drop/ml

- Spike protector, 40mm, PE
- Spike, 20 drop/ml, ABS
- Vented air cap VAC02, PVC, Hydrophobic membrane 0.5µm
- Chamber, Φ15.5x64mm, PVC
- Liquid lock membrane, 8µm

N7F

64无滤网滴瓶头部, 20滴
64mm, Non-vented assembled chamber, 20 drop/ml

- Spike protector, 40mm, PE
- Spike, 20 drop/ml, ABS
- Chamber, Φ15.5 x 64mm, PVC.
- Liquid lock membrane, 8µm

成品输液器套装--止液型

We Can Customize Liquid Lock IV Sets For You



Gravity IV Set

- Standard IV set for infusion of drugs and liquids
- Liquid lock drop chamber, 20 drops/ml, filter 8µm
 - Roller clamp
 - Needle free Y injection site
 - Male luer lock




Pump IV Set

- Pump IV set for infusion of drugs and liquids
- Liquid lock drop chamber, 20 drops/ml, filter 8µm
 - Roller clamp
 - Y-port injection site, latex-free
 - Slide clamp
 - Male luer lock



PVC Free IV Set

- IV set for infusion of drugs and liquids
- Liquid lock drop chamber, 20 drops/ml, filter 8µm
 - Roller clamp
 - Y-port injection site, latex-free
 - Male luer lock
 - PVC free



Light Resistant IV Set

- Light Resistant IV set for infusion of drugs and liquids
- Liquid lock drop chamber, 20 drops/ml, filter 8µm
 - Roller clamp
 - Flashback
 - Male luer lock

DEHP Free
PVC Free