




Angiographic Accessory Products 2021



Coronary angiography is an important means for cardiologists to diagnose coronary heart disease and interventional stent.


In clinical imaging, a catheter is used to run along the invasive entrance of the radial or femoral artery to the coronary artery opening of the heart. The contrast medium is then injected into the coronary artery so that the internal shape of the coronary artery can be shown on the X-ray to facilitate the vascular surgeon's evaluating of the disease.

Infusion sets is an important medium for the circulation of contrast agent and normal saline during the surgery of angiography.

Phone:(86)21-5774 3953
www.bqplusmedical.com

BQ+ Medical Co., Ltd
上海必趣医疗科技有限公司

Angiographic Accessory Products



Infusion set - Contrast Administration set (with IV burette chamber)




- * Contrast medium dosage can be controlled
- * Various sets are optional
- * Customized
- * DEHP FREE/ LATEX FREE

Item	Specifications	Notes
IV burette sets (without drip chamber, 180cm) Tube ID3.0mm*OD4.0mm	individually packed with peel pouch 10pcs/case, 100pcs/ctn	layout customized EO sterilized
IV drip burette sets (with drip chamber, 180cm) Tube ID3.0mm*OD4.0mm	individually packed with peel pouch 10pcs/case, 100pcs/ctn	layout customized EO sterilized

Angiographic Accessory Products


Infusion set - Contrast Administration set



- * Various sets are optional
- * Customized
- * DEHP FREE/ LATEX FREE

Item	Specifications	Notes
IV sets with spike (Tube ID3.0mm*OD4.0mm,180cm)	individually packed with peel pouch 25pcs/case, 300pcs/ctn	layout customized EO sterilized
IV sets with spike (Tube ID3.9mm*OD5.5mm,180cm)	individually packed with peel pouch 25pcs/case, 300pcs/ctn	layout customized EO sterilized
IV sets with drip chamber (vented,Tube ,180cm) Tube ID3.0mm*OD4.0mm	individually packed with peel pouch 25pcs/case, 300pcs/ctn	layout customized EO sterilized
IV sets with drip chamber (non-vented,180cm) Tube ID3.0mm*OD4.0mm	individually packed with peel pouch 25pcs/case, 300pcs/ctn	layout customized EO sterilized

Infusion set - Fluid Administration Set



- * Various sets are optional
- * Customized
- * DEHP FREE/ LATEX FREE

Item	Specifications	Notes
IV sets with drip chamber (vented、180cm) Tube ID3.0mm*OD4.0mm	individually packed with peel pouch 25pcs/case, 300pcs/ctn	layout customized EO sterilized
IV sets with drip chamber (non-vented、180cm) Tube ID3.0mm*OD4.0mm	individually packed with peel pouch 25pcs/case, 300pcs/ctn	layout customized EO sterilized