

Injection Site 加药件



ST103-Y

嵌入式三通加药口合成胶
Y-injection site, latex-free

- Size: length=35mm, with easy handle
- Material: ABS, Isoprene rubber



ST103-OR

嵌入式三通加药口合成胶
Y-injection site, latex-free

- Size: length=35mm, with easy handle
- Material: ABS, Isoprene rubber



ST104-Y

嵌入式三通加药口合成胶
Y-injection site, latex-free

- Size: length=32.4mm, with easy handle
- Material: PC, Isoprene rubber
- Pre-pierced available



ST102-IR

嵌入式三通加药口合成胶
Y-injection site, latex-free

- Size: length=35mm
- Material: ABS, Isoprene rubber



ST108-IR

嵌入式三通加药口合成胶
Y-injection site, latex-free

- Size: length=38.9mm
- Material: PC, Isoprene rubber
- Pre-pierced available



ST113-IR

嵌入式三通加药口合成胶
Y-injection site, latex-free

- Size: length= 35mm
- Material: ABS, Isoprene rubber

STY08-IR

透明肝素帽
Transparent male luer lock injections site

- Heparin cap, latex-free
- Material: ABS, Isoprene rubber



STY09-IR

八角透明肝素帽
Transparent male luer lock injections site

- Heparin cap, latex-free
- Material: ABS, Isoprene rubber



STY07-IR

合成胶肝素帽
Male luer lock injection site

- Heparin cap, latex-free
- Material: ABS, Isoprene rubber



STY10-IR

透明肝素帽
Transparent male luer lock injections site

- Heparin cap, latex-free
- Material: PC, Isoprene rubber
- Pre-pierced available



STY11-IR

合成胶肝素帽
Male luer lock injection site

- Heparin cap, latex-free
- Material: PC, Isoprene rubber
- Pre-pierced available



AC02

乳胶泡合成胶
Flashback, latex-free

- Size: length=45mm
- Material: Isoprene rubber



Roller Clamp 滚轮流量调节器

RC101-B

长尾调节器
Precision V Roller clamp

- Size: length=56mm
- For tubing O.D.4.0mm
- Material: ABS



RC101-B

RC102-B

短尾调节器
Precision V Roller clamp

- Size: length= 56mm
- For tubing O.D.4.0mm
- Material: ABS



RC102-B

RC103

小调节器
Precision V Roller clamp

- Size: length= 45mm
- For tubing O.D.4.0mm
- Material: ABS



RC102-C

短尾调节器
Precision V Roller clamp

- Size: length=56mm
- For tubing OD 4.0-4.2mm
- Material: ABS



RC102-C

RC108

止液调节器
Roller clamp used on air stop infusion set

- Size: length= 56mm
- For tubing O.D.4.0mm
- Material: ABS

