

Safety ENFit Enteral Nutrition Feeding Set

喂食系列


Enteral Feeding Bags & Enteral Feeding Sets 喂食系列



Enteral Feeding Bags

喂食袋套装



- ① Feeding bag(Gravity/Pump, 500/1000/1200ml)
- ② PVC tubing (DEHP/Non-DEHP)
- ③ Roller clamp
- ④ IV chamber
- ⑤ Silicon tubing(Optional)
- ⑥ Silicon connector(Optional)
- ⑦ Feeding Y connector with cap
- ⑧ Enfit lock female luer with cap
- ⑨ Caution tag
- ⑩ Patient label



Enteral Feeding Sets

喂食管套装


- ① Spike with cap
- ② PVC tubing (DEHP/Non-DEHP)
- ③ Roller clamp
- ④ IV chamber
- ⑤ Silicon tubing(Optional)
- ⑥ Silicon connector(Optional)
- ⑦ Feeding Y connector with cap
- ⑧ Enfit lock female luer with cap
- ⑨ Caution tag
- ⑩ Patient label



ENFit Extension Sets


喂食延长管

- ① Enfit lock female luer with cap
- ② PVC tubing (DEHP/Non-DEHP)
- ③ Enfit lock male luer with cap
- ④ Tubing OD 3.0/4.0mm



Safety ENFit Syringes & ENFit Components

喂食注射器&喂食配件



ENFit Syringes

喂食注射器

Volume	Graduation	Applications	Optional Package
1ml, 3ml, 5ml, 10ml	0.2ml	Medicine, Flush	Blister Pack
20ml, 35ml, 60ml	1.0ml	Hydrate, bolus feed	Back Pack

ENFit Components 喂食配件



EN-CP02

Spike with protector



EN-YS01

Y-connector with cap

- Y-connector EN-ST115Y
- Adapter ring ENS01
- Female cap EN-HT132



EN-CN07

Adaptor for syringe



RC101-B

Roller clamp



EN-MH01

ENFit lock male connector with cap

- Male Adapter EN-MLL15
- Adapter ring ENS01
- Female cap EN-HT132
- For tubing OD 3.0/4.0mm



EN-FH01

ENFit lock female connector with cap

- Female Connector EN-FLL05
- Adapter ring ENS01
- Male cap EN-HT133
- For tubing OD 3.0/4.0mm



EN-CP05

ENFit transition connector



EN-CP04

Christmas tree connector