

Luer Lock 鲁尔接头



FLL01

母螺旋
Female Luer Lock

- Size: For tube OD 2.5/3.0/4.0mm
- Material: ABS/PC
- Cap: HT105, PE



FLL02

母螺旋
Female Luer Lock

- Size: For tube O.D. 3.0/4.0/5.0/5.4mm
- Material: PC/ABS



FLL05

母螺旋
Female luer lock

- Size: For tubing OD 4.0mm
- Material: PC



FLL04

母螺旋
Female luer lock

- Size: For Tubing I.D. 4.0/4.7/6.4mm
- Material: PC



MLL02

螺旋座
Male luer lock

- Size: length=32mm, for flashback
- Material: ABS
- Cap: HT102, PE



MLL03

内外扣管式螺旋座
Male luer lock

- Size: length=40mm, for tube 3x4.1mm
- Material: ABS
- Cap: HT103, PE



MLL04

公螺旋
Male Luer Lock for extension

- Size: For tube OD 2.5/3.0/4.0mm
- Material: ABS/PC
- Cap: HT104, PE



MLL10

公螺旋
Male Luer Lock

- Size: For tubing I.D. 4.0/4.7/6.4mm
- Material: PC



MLL13

公螺旋
Male luer lock

- Size: For Tubing O.D. 4.0mm
- Material: ABS/PC



DMC-01

混药接头
Fluid Dispensing Connector

- Female to female ends
- Sterile, blister package



HT110-2N

护帽
Cap

- Material: PE/PP
- L= 19mm
- Non-vented



CN01/CN13

三通
Connectors

- Size: For tube O.D. 2.3/4.0mm
- Material: PVC

Luer Lock 鲁尔接头



MLL05

三联件/可脱卸螺旋座
Detachable spin luer lock

- Size: for tube OD 4.0mm
- Material: ABS
- Cap: HT106, PE



MLL05-OR

避光三联件/可脱卸螺旋座
UV Detachable spin luer lock

- Size: for tube OD 4.0mm
- Material: ABS
- Cap: HT106, PE



MLL05-B

三联件/可脱卸螺旋座
Distachable spin luer lock

- Size: for tube OD 4.0mm
- Material: ABS
- Cap: HT130, PP, blue/green/red
With hydrophobic membrane 0.2/1.2/3.0µm



MLL12

三联件/可脱卸螺旋座
Detachable spin luer lock

- Size: for tubing OD 2.0/2.2/2.4/2.5mm
- Material: ABS/PC
- Cap: HT107, PE



MLL14

不可拆卸/可旋转螺旋座
Undetachable spin luer lock

- Size: for tubing OD 3.6/4.0mm
- Material: PC
- Cap: HT106



MLL104

螺旋座
Male Luer Lock

- Material: ABS
- For Tube OD 2.9/4.0MM



MLS01

带翅膀三通
Male luer slip

- Size: length=30mm, for flashback
- Cap: HT107, PE



MLS02

小三通
Male Luer slip

- L=29mm
- Size: For Flashback
- Cap: HT107, PE



MLS03

三联件三通
Male Luer slip

- L=28mm
- Size: For tube OD 3.0/4.0 mm
- Cap: HT107, PE

Cap 护帽



HT151

护帽
Cap

- Material: ABS
- For stopcock



HT152

护帽
Cap

- Material: PE
- For stopcock



HT153

带孔护帽
Cap With Hole

- Material: PP