

Spike & Back Check Valve 塑针&单向阀

Spike 塑针



ST03

双头塑针
Spike

- Double Ended Transfer Spike
- Material: ABS & PE



ST02

塑针
Spike

- For tubing O.D.4.1mm
- Material: ABS

ST18

塑针
Vented Spike

- For tubing O.D.4.0mm
- Material: ABS



VST01

塑针
Vented Spike

- For tubing O.D.4.0mm
- Material: ABS



Back Check Valve 单向阀

BCV


单向阀
Back Check Valve (One-Way Valve)

- Back Pressure: > 300kPa
- Flow rate: 100ml/min
- Latex Free



Ret.No.	Product Description
BCV06	Female-Male Ends
BCV08	Male-Tube Ends 2.3*3.5mm (I.D.xO.D.)
BCV10	Female-Tube Ends 2.3*3.5mm (I.D.xO.D.)
BCV12	Tube-Tube Ends 2.3*3.5mm (I.D.xO.D.)

High Pressure Manifold & Stopcock 高压联板&高压三通阀



High Pressure Manifold

高压联板


- No Leakage under 500psi
- Smooth rotation, easy use
- Different color for choice

	Ref.No.	Product Description
2 Ports	HPM201	"ON" type with rotating adapter, Right port orientation
	HPM202	"OFF" type with rotating adapter, Right port orientation
3 Ports	HPM301	"ON" type with rotating adapter, Right port orientation
	HPM302	"OFF" type with rotating adapter, Right port orientation
	HPM303	"ON" type with pressure line, Right port orientation
	HPM307	"ON" type with rotating adapter, left port orientation
	HPM308	"OFF" type with rotating adapter, left port orientation
4 Ports	HPM401	"ON" type with rotating adapter, Right port orientation
	HPM402	"OFF" type with rotating adapter, Right port orientation

High pressure Stopcock

高压三通阀

- No Leakage under 1000psi Pressure
- Different color for choice
- Smooth rotation, easy use



Ref.No.	Product Description
HPSK0106	"ON" Type, Rotating Adapter, No Leakage under 1000psi Pressure
HPSK0108	"ON" Type, Rotating Adapter, No Leakage under 500psi Pressure
HPSK0109	"OFF" Type, Rotating Adapter, No Leakage under 1000psi Pressure
HPSK0110	"OFF" Type, Rotating Adapter, No Leakage under 500psi Pressure