



Blood Transfusion Set 输血器套装




Blood Set
普通输血器



Air-vented Blood Set

带空气塞输血器



High Pressure Extension Set 高压延长管

High Pressure Tubing

高压延长管

- Multiple Length And Connector Types For Choice
- PU Tubing With Nylon Liner For 1200psi
- PVC Tubing For 500psi




Female-Rotaing Adapter



Female-Male Luer Lock



3-Way Stopcock-Male Luer Lock



Ret.No.	Type	Product Description	Ret.No.	Type	Product Description
HPEXT101	30cm	Female-Male Luer Lock No Leakage Under 500psi Pressure	HPEXT501	30cm	Rotaing Adapter No Leakage Under 1200psi Pressure
HPEXT102	60cm		HPEXT502	60cm	
HPEXT103	90cm		HPEXT503	90cm	
HPEXT104	120cm		HPEXT504	120cm	
HPEXT105	150cm		HPEXT505	150cm	
HPEXT106	180cm				
HPEXT301	30cm	Female-Male Luer Lock No Leakage Under 1200psi Pressure	HPEXT201	30cm	3-Way Stopcock-Male Luer Lock No Leakage Under 40psi Pressure
HPEXT302	60cm		HPEXT202	60cm	
HPEXT303	90cm		HPEXT203	90cm	
HPEXT304	120cm		HPEXT204	120cm	