



Invasive Blood Pressure Transducer (IBP) 2021




In the treatment of critically ill patients process in recent years (for example: shock patients, cardiac surgery, and organ transplantation, etc.), it is necessary to maintain several pressure parameters such as arterial pressure and venous pressure in human body, for the continuous dynamic monitoring, invasive pressure sensors are widely used in clinical practice.

Phone:(86)21-5774 3953
www.bqplusmedical.com

BQ PLUS MEDICAL CO., LTD

Invasive Blood Pressure Transducer (IBP)



BQ Products




Advantages:

Plastic needle puncture force is enhanced, good toughness, better bonding with PVC chamber, no leakage

Invasive Blood Pressure Transducer (IBP)

Infusion set--single way





* Department of use: operating room, anesthesiology, all kinds of ICU, catheter room

* Widely used in China

Configuration	chamber+roller clamp+tubing+male luer
Tubing length	1500mm
Chamber length	64mm/55mm
Package	1 set/, 25sets/PE bag, 450/ctn

Infusion set--double way




* Department of use: operating room, anesthesiology, all kinds of ICU, catheter room

* Blood pressure can be measured at two sites simultaneously. Widely used abroad

Configuration	chamber+roller clamp+connector+tubing+male luer
Tubing length	1730±50mm/350±30mm
Chamber length	64mm/55mm
Package	1 set/, 15/PE bag, 300/ctn

Invasive Blood Pressure Transducer (IBP)

Infusion set -- triple way



- * Department of use: operating room, anesthesiology, all kinds of ICU, catheter room
- * Blood pressure can be measured at three sites simultaneously. Widely used abroad

Configuration	chamber+roller clamp+connector+tubing+male luer
Tubing length	1730±50mm/350±30mm
Chamber length	64mm/55mm
Package	1 set/, 10/PE bag, 200/ctn

Non-compressible tubing



- * Under 15N static tension 15S does not fall off or disconnect

Configuration	female luer+tubing+male luer
Long tube length	1200mm / 1500mm
Short tube length	300mm
Pressure	500psi

Invasive Blood Pressure Transducer (IBP)

Stopcock (one-way & three-way)



- * Pressure: 300kpa water pressure
- * handle is stopped every 45 degrees angle
- * luer lock rotates 360° to prevent tubing kinking

- * Single handle (OFF), easy to operate for nurses widely used in the market



The handle of the stopcock is smooth , easy to operate, which is helpful to discharge liquid and blood sample collection in the process of monitoring

Accessories



- * With all accessories(such as roller clamp, robert clamp,pinch clamp,caps with or without holes,connectors
- * Various colors
- * Customized