




Micro Pump Extension Tube

2021




The micro pump extension tube can also be used as an extension tube for the micro syringe pump. It is used in conjunction with the micro syringe pump. Pediatric patients need to use drugs and fixed or fixed speed drugs. After assembling a syringe and an intravenous infusion needle, it can be supplied to the human body in the form of infusion. It can be used for liquid injection, and can also be used to lengthen the pipeline when the length of the pipeline is not enough. The product is for one-time use.

The micro injection pump can accurately, uniformly and continuously inject drugs into the human body, reasonably adjust the injection speed of the drug, strictly control the drug dosage, and ensure the best effective concentration of the drug. It is used for continuous infusion of various urgently needed drugs, such as All kinds of vascular activity drugs, anti-arrhythmic drugs, electrolyte solutions, anesthetics, pediatric drugs infusion, etc., to reduce the occurrence of complications, and its simple operation can reduce workload when rescuing critically ill patients, improve work efficiency and accuracy , Safely and effectively cooperate with doctors for rescue.

Phone:(86)21-5774-3953
www.bqplusmedical.com

BQ+ Medical Co., Ltd


Micro Pump Extension Tube



The micro pump extension tube provided by our company is mainly divided into : normal type, pressure type, light-resistance type and extension sets.

1. The complete type of ordinary extension tube is composed of luer lock, pipeline, clamp, presicion filter and protective cap.
2. The raw materials used are ABS, PE and PVC.
3. Performance: sterile, pyrogen-free, the pore size of the liquid medicine filter is 0.2/1.2/5µm.

Pressure Type



Reorder Code	Description	Features	Approx. Overall Length	Units/ Case
YV-1	Micro Bore Extension with Slide Clamp, Male&Female Luer Lock.	Non-DEHP Formulation, Latex-Free, withstand 250Kpa pressure for 10 seconds, and there should be no gas leakage.	≥50cm	50


Normal Type







Reorder Code	Description	Features	Approx. Overall Length	Units/ Case
YV-2	Micro Bore Extension with back check valve, Slide Clamp, Male&Female Luer Lock.	Non-DEHP Formulation, Latex-Free, withstand 200Kpa pressure for 10 seconds, and there should be no gas leakage.	≥50cm	50

PERIPHERAL EXTENSION SETS

Precision Filter Extension Sets








BQ+ Code	Product Description	Size	Units/case
YV-3	Micro Bore Extension with Micron Filter, slide Clamp, Male&Female Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	25
			
YV-4	Micro Bore Extension with Baby droplet Micron Filter, slide Clamp, Male&Female Luer Lock. Non-DEHP Formulation, Latex-Free, 1.5*3.0mm tubing IN/OUT	≥50cm	25
			
YV-5	Micro Bore Extension with Micron Filter, slide Clamp, needle free valve, Male & Female Luer Lock. Non-DEHP Formulation, Latex-Free, PP, ABS and Silicone.	≥50cm	25
			
YV-6	Micro Bore Extension with Micron Filter, robert Clamp, injection site, Male&Female Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	25
			

PERIPHERAL EXTENSION SETS

Needle Free Extension Sets







BQ+ Code	Product Description	Size	Units/case
YV-7	Micro Bore Extension with needle free valve,slide clamp, Male&Female Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-8	Micro Bore Extension with two needle free valves, two slide clamps,Bifurcated connector, Male & Female Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-9	Micro Bore Extension with needle free valve, female luer lock, two slide clamps,Bifurcated connector, Male&Female Luer Lock Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-10	Micro Bore Extension with luer lock,Y-needle free valve, slide clamps,Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-11	Micro Bore Extension with luer lock, two Y-needle free valves,slide clamps,Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			

PERIPHERAL EXTENSION SETS

Stopcock Extension Sets







BQ+ Code	Product Description	Size	Units/case
YV-12	Micro Bore Extension with Stopcock, Luer Lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-13	Micro Bore Extension with Stopcock, Luer Lock,slide clamp, luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-14	Micro Bore Extension with Stopcock, Y Needle free valve, Luer Lock,slide clamp, luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-15	Micro Bore Extension with two Stopcocks,luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			

PERIPHERAL EXTENSION SETS

Back Check Valve Extension Sets



BQ+ Code	Product Description	Size	Units/case
YV-16	Micro Bore Extension with two robert clamps, back check valves, Bifurcated connector, luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-17	Micro Bore Extension with two back check valves, Bifurcated connector, luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-18	Micro Bore Extension with two back check valves, needle free valve, luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			
YV-19	Micro Bore Extension with slide clamp, two back check valves, two needle free valves, luer lock. Non-DEHP Formulation, Latex-Free	≥50cm	50
			

PERIPHERAL EXTENSION SETS

Protected Spiral Extension Sets



BQ+ Code	Product Description	Size	Units/case
YV-20	Protected Spiral Extension Line, 3-layer extension tube, back check valve, Female / Male Luer-lock. Non-DEHP Formulation, Latex-Free,	≥200cm	30

Address: No.18,Cheye Road,Songjiang District,
Shanghai 201611 China.

Email: info@bq-medical.com

Tel: (86)21-5774 3953

Web: www.bqplusmedical.com