


BQ+Needle Free Valve

无针加药件



Advantages:

- MRI compatibility
- Blood compatibility
- Smooth and tight seal, easy to disinfect thoroughly
- Good leakproofness
- Good self-resilience

Item	Functional Activation	Priming Volume	Flow Rate	Leakproofness	Material	Dimension(mm)	Weight(g)
NFV-L16	Over 200 times	0.06ml	2000ml/min	No leakage under 200kpa by air pressure & hydraulic pressure	PP, ABS and Silicone	35*11*11	1.7
NFV-S16		0.09ml	850ml/min			37*11*11	1.7
NFV-Y16		0.47ml	850ml/min			48*11*11	2.3

Finished IV Set with Needle Free Valve

成品输液套装--带无针加药件



Needle Free Valve

Blister Package
OEM available



Extension Set

Extension set for infusion and tranfusion
·Female luer lock
·Needle free valve
·Male luer lock
·Tubing Length: 15/30/60/90/120cm



Gravity IV Set

Standard IV set for infusion of drugs and liquids
·Drop chamber, 20 drops/ml, filter 15um
·Roller clamp
·Needle free valve
·Male luer lock



IV Burette Set

IV set for Burette set
·Vented spike
·Roller clamp
·100/150ml Burette chamber(optional)
·Roller clamp
·Flow regulator
·Needle free valve
·Male luer lock

