


BQ+ High Flow Precision Filter

大流量精密过滤器



- Features:**
- No leakage under 300Kpa 15S
 - High flow / fast priming available
 - Remove bubble
- Material:**
- ABS body
 - PES membrane
 - PTFE membrane



Description	Color	Aperture	Filter efficiency	Minimum flow reate	Inlet & outlet	Effective filtration area
Adult Square filter High flow sery Air-Vented	Clear /Amber	0.2um	Sterilizing-grade	>320mL/10min	2.9*4.05mm Tubing IN/OUT FLL IN/RMLL OUT	12 cm ²
		1.2um	≥90%	>1000mL/10min		
		5um	≥90%	>1000mL/10min		
Kid Square filter Air-Vented	Clear /Amber	0.2um	Sterilizing-grade	>150mL/10min	2.9*4.05mm Tubing IN/OUT FLL IN/RMLL OUT	6 cm ²
		1.2um	≥90%	>900mL/10min		
		5um	≥90%	>1000mL/10min		
Baby Droplet filter Air-Vented	Clear /Amber	0.2um	Sterilizing-grade	>40mL/10min	1.5*3.0mm Tubing IN/OUT FLL IN/RMLL OUT	1.5 cm ²
		1.2um	≥90%	>300mL/10min		
		5um	≥90%	>550mL/10min		

Finished IV Set with High Flow Precision Filter

成品输液套装—带大流量精密过滤器



Extension Set

- Extension set for infusion of drugs and liquids
- Female luer lock
 - Precision filter 0.22/1.2/5um(optional)
 - Male luer lock
 - Tubing Length: 15/30/60/90/120cm



Gravity IV Set

- Standard IV set for infusion of drugs and liquids
- Drop chamber, 20 drops/ml, filter 15um
 - Roller clamp
 - Precision filter 0.22/1.2/5um(optional)
 - Male luer lock



Pump IV Set

- Pump IV set for infusion of drugs and liquids
- Drop chamber, 20 drops/ml, filter 15um
 - Roller clamp
 - Precision filter 0.22/1.2/5um(optional)
 - Needle free valve
 - Male luer lock



IV Burette Set

- IV set for Burette set
- Vented spike
 - Roller clamp
 - 100/150ml Burette chamber(optional)
 - Roller clamp
 - Precision filter 0.22/1.2/5um(optional)
 - Needle free valve
 - Male luer lock

